

## DAILY CLEANING AND CARE OF GEL LINER



Remove liner from limb and turn inside out

With the GEL SIDE OUT, using a mild soap (i.e. ivory soap) thoroughly clean the surface of your liner. Then RINSE ALL the soap and residue from the liner with water.



Pat the liner dry with a towel.



Place the liner FABRIC SIDE OUT on the provided drying stand. Then clean your limb with soap and water and rinse well.

Never leave the liner with the gel side out



**Proper hygiene of your ALPHA liner and Limb are very important and both should be cleaned at the end of each day of wear.**